

# **Simple Annotation of Metatranscriptomes through Sequence Annotation (SAMSA): A Pipeline for Metatranscriptome Analysis**

SAMSA version 1.0.0

Created by Sam Westreich

[stwestreich@ucdavis.edu](mailto:stwestreich@ucdavis.edu)

<http://www.github.com/transcript/>

## **Table of Contents**

Introduction to SAMSA (and core dependencies)	3
A note for Windows users	4
SAMSA Pipeline Use: Streamlined wrapper script for pre-annotation	5
SAMSA Pipeline Use: Detailed pre-annotation steps	7
SAMSA Pipeline Use: Detailed post-annotation steps	10
Example files and results	17
Frequently Asked Questions (FAQs)	22
Contact information	24

## Introduction to SAMSA (and core dependencies)

Hello, and thank you for considering using SAMSA! This pipeline is designed to take raw sequencing reads from an Illumina run and match those reads against a reference database, and then examine the annotations for significant changes in either organism activity, or in the expression levels of various transcripts.

Although metatranscriptomics—the study of diverse microbial population activity based on RNA-seq data—is rapidly growing in popularity, there are limited options for biologists to analyze this type of data. Current approaches for processing raw metatranscriptome data rely either on restricted databases, a personal analysis server, or use metagenome-based approaches that have not been fully evaluated for use in processing metatranscriptomic datasets. I have created a new bioinformatics pipeline, SAMSA, designed specifically for metatranscriptome dataset analysis, which runs either in-house or in conjunction with Metagenome-RAST (MG-RAST) servers. Designed for use by researchers with relatively little bioinformatics experience, SAMSA offers a breakdown of metatranscriptome activity by organism or transcript function, and is fully open source.

SAMSA was constructed in a Mac OS X environment, and is designed to run in Unix systems. The majority of the pipeline is written in Python, with R scripts used in the final steps for performing statistical comparisons and generating visual graphs of the metatranscriptome composition.

### Version

This is SAMSA version 1.0.0, created in winter 2016.

### Dependencies

- Unix/Mac OS environment. If running on a Windows computer, Cygwin (<https://www.cygwin.com/>) may be used to simulate a Unix command line, although all files must be moved into the Cygwin directory.
- Python. This pipeline was constructed to run on Python 2.7.
- Python modules:
  - sys
  - os
  - subprocess
  - glob
  - time
  - gzip
  - operator
  - commands
- R, or more preferably, RStudio. Analysis R scripts use the following R packages, which can ALL be installed by running the included script `package_install.Rmd` in R:
  - DESeq2

- ggplot2
- gridExtra
- scales
- knitr
- reshape
- Plyr
- Trimmomatic, a flexible read-trimming tool for NGS data. Documentation and download links for Trimmomatic may be found here: <http://www.usadellab.org/cms/?page=trimmomatic> .
- FLASH, Fast Length Adjustment of SHort reads, a tool for merging paired end reads. If the sequencing data is not paired end, this tool is not necessary - although paired end data is more accurate. FLASH download and documentation links may be found here: <https://ccb.jhu.edu/software/FLASH/> .
- A working internet connection. SAMSA needs to both upload and download files from MG-RAST's servers, which may require several hours of uninterrupted access.
- Memory: all SAMSA commands have been shown to work on a laptop with only 8 Gb of RAM, but more ram will increase processing speed.
- An MG-RAST account. MG-RAST can be accessed at <http://metagenomics.anl.gov/>, and creating an account is free. SAMSA was built to interface with MG-RAST version 3.6, which is current as of May 2016.

## A note for Windows users

SAMSA is designed to run in a Mac OS/Unix environment. However, don't despair! With a couple tweaks, SAMSA can still run on a Windows machine.

First, for all command line steps, Cygwin (<https://www.cygwin.com/>) is necessary. This program emulates a Unix environment. Simply install, and move all SAMSA programs into the Cygwin directory (generally labeled as home/). Cygwin comes with Python installed. Both FLASH and Trimmomatic will need to be moved into the Cygwin home directory as well.

For the R programs, RStudio is available for Windows (<https://www.rstudio.com/products/rstudio/download/>), and will allow for all R commands to be carried out.

In addition, running all lines in the script labeled package\_install.Rmd will install all needed R packages for the pipeline.

## Output files

Note that, at each stage in the pipeline, SAMSA generates output files, which serve as the inputs for the next stage in data processing. This allows for the pipeline to be paused at any step, and different steps may be run on different machines by transferring the files and pipeline programs.

Now, ready to get started?

## SAMSA Pipeline Use: Streamlined wrapper script for pre-annotation

For ease of use, several wrapper programs exist to run the various steps in this pipeline with a minimum of interacting with the command line. These wrapper scripts are entitled:

1. SAMSA\_pre\_annotation\_pipeline.py: this wrapper script performs the following steps:
  - a. Read cleaning and adaptor removal with Trimmomatic
  - b. Paired-end file aligning with FLASH (if the “paired” option is enabled)
  - c. Uploading of files to MG-RAST
  - d. Sequence statistics computation with MG-RAST
2. SAMSA\_post\_annotation\_pipeline.py: this wrapper script performs the following steps:
  - a. Downloading organism and functional annotations from MG-RAST
  - b. Running SAMSA analysis scripts on downloaded annotation files
  - c. Cleans up summary output files for import into R for analysis

### Pre-annotation wrapper program

SAMSA\_pre\_annotation\_pipeline.py has the following usage options:

#### REQUIRED:

- Ends (1 or 2) Specifies whether input sequence files are single or paired end
- A (string) Authorization key generated by MG-RAST
- D (path) Path to folder containing raw sequence files to be processed

#### OPTIONAL:

- Q Enables ‘quiet mode’; no output printed to the console.
- F (path) Path to FLASH app, if not installed in /Applications/FLASH-1.2.11/flash\* .
- T (path) Path to Trimmomatic.jar, if not in

/Applications/trimmomatic-0.33/trimmomatic.jar .

Example pre-annotation command:

```
$ python SAMSA_pre_annotation_pipeline.py -A MG-RAST-KEY -D  
~/path/to/sequences/ -Ends 2 -F ~/path/to/FLASH/program* -T ~/path/to/trimmomatic.jar
```

### Submitting files on MG-RAST for annotation

After the SAMSA\_pre\_annotation\_pipeline.py wrapper has finished, the files are trimmed, aligned (if paired end), uploaded to MG-RAST’s inbox, and MG-RAST has computed sequence stats for each file.

To submit these files, log in to MG-RAST and select “Upload”, and then “Next” at the bottom of the page. Select all files to be submitted, assign a project name, and make sure to DESELECT “dereplication.” After the files have been submitted, their progress can be viewed by selecting “Browse” at the top of the page, and then “In Progress” on the left side.

## Downloading annotations from MG-RAST

The best approach for downloading annotation files from MG-RAST is to use the MG-RAST\_API\_downloader.py, part of the SAMSA tools. This python script has the following necessary command line arguments:

### REQUIRED:

- S (string)      Annotation source: RefSeq, UniProt, KEGG, etc.
- D (string)      Data type: organism, function, or ontology
- A (string)      MG-RAST authorization key
- I (number)      Annotation ID number; found on the MG-RAST results page
- O (name)      Save name for this annotation file (will be created in current directory).

### OPTIONAL:

- Q              Enables quiet mode.
- usage          Prints these usage options to the command line output and then exits.

These commands can also be seen by running:

```
$ python MG-RAST_API_downloader.py -usage
```

A sample MG-RAST\_API\_downloader.py command:

```
$ python MG-RAST_API_downloader.py -S RefSeq -D organism -I 459902.3 -A  
UyTjAgVxLPwwdGDj73xhJ9w9C -O output_metatranscriptome.tab
```

Because this download program can be run from a single line of code, it is possible to create a very simple bash script containing download program commands for each of the metatranscriptomes to be downloaded. This allows for a “queue” of downloads.

NOTE: Given the size of these annotation files, the download process can take several hours. It's advisable to run this process in a separate shell using the screen command (<https://www.gnu.org/software/screen/manual/screen.html>), allowing for logging off while the download continues to run.

**Once these annotation files are downloaded, proceed to Step 5, in later instructions, for analysis.**

## **SAMSA Pipeline Use: Detailed pre-annotation steps**

Steps 1-3 are performed automatically if using the SAMSA\_pre\_annotation\_pipeline.py, as described above. However, they can each be performed manually based off of the steps described here.

### **Step 1: Preprocessing**

Raw files must be processed using Trimmomatic to remove adaptor contamination and low-quality reads before they can be submitted to MG-RAST.

For details using Trimmomatic, refer to the manual:

[http://www.usadellab.org/cms/uploads/supplementary/Trimmomatic/TrimmomaticManual\\_V0.32.pdf](http://www.usadellab.org/cms/uploads/supplementary/Trimmomatic/TrimmomaticManual_V0.32.pdf)

*If using paired-end reads:*

For each sample, two files should be present, often designated as R1 and R2. The R1 file contains all forward reads for that sample, while the R2 file contains all reverse reads.

FLASH (<http://ccb.jhu.edu/software/FLASH/>) is recommended for joining these two files together to create a single file of all combined reads. FLASH can be downloaded through SourceForge. Once downloaded, run the following commands:

```
$ tar xvf FLASH-1.2.11.tar.gz
$ cd FLASH-1.2.11/
$ make
```

FLASH is run from the command line, using the following command:

```
$ ./flash MATES_1.fastq MATES_2.fastq
```

To see FLASH options, use the following command:

```
$ ./flash --help | less
```

### **Step 2: Generating an authorization key and uploading to MG-RAST**

Once files have been cleaned (and joined if paired-end), the next step is to upload the file to MG-RAST for annotation.

If you do not have an MG-RAST, it is necessary to make an account at

<http://metagenomics.anl.gov/?page=Register>. If you do have an account, you will need a current authorization key (generated within the last ten days). The authorization key can be found at <http://metagenomics.anl.gov/?page=AccountManagement> by selecting “Preferences”.



If you generate a new authorization key, be sure to save the key by selecting “Set Preferences” at the bottom of the page.

**MG-RAST**  
metagenomics analysis server

Sam Westreich

### Manage Preferences

#### MANAGE PREFERENCES

Your preferences are divided into categories, which generally represent pages. If you have not yet chosen a preference for a certain setting, the default value will be used. You can come to this page and change your preferences at any time.

##### BROWSER

do not display incorrect browser popup

set preferences

##### FUNDING

funding sources

set preferences

##### WEB SERVICES

authentication key for web services

webkey termination date

generate new key

Note: Creating a new key and clicking 'set preferences' will render the previous key deprecated. Your key will be valid for a limited time only (see webkey termination date). You can generate a new key with a new termination date at any time.

set preferences

Although files can be uploaded to MG-RAST through the point-and-click interface, it is simpler and more rapid to use the API interface. To upload files through the MG-RAST API:

1. Navigate to the folder on the Unix system where the files are located.
2. Run `upload_MG-RAST.py` (part of the SAMSA tools). The following options must be specified:
  - a. `-A authorization_key`
  - b. `-F file_to_be_uploaded`
  - c. For further usage options, run: `upload_MG-RAST.py -usage`
3. Files may be uploaded in zipped form (.gz) for faster transfer.

Example `upload_MG-RAST.py` command:

```
$ python upload_MG-RAST.py -A eZYmtTzx5yrYn8pDQx4kUFxC4 -F sequences.fastq
```

### Step 3: Submitting to MG-RAST

Once all files are uploaded to MG-RAST's inbox, access the inbox at <http://metagenomics.anl.gov/Html/mgmainv3.html?mgpage=upload>. MG-RAST must compute sequence stats on each file before it can be submitted. If files were uploaded in zipped format, they must be unzipped in the inbox.

To submit files for processing, click the “Next” button at the bottom right. When submitting a metatranscriptome, be sure to ensure that under “Pipeline Options”, “Dereplication” is UNCHECKED (disabled). *Warning: failing to disable dereplication will result in the removal of identical sequences and the loss of expression data!*

#### 4. choose pipeline options

**assembled** ☐ Select this option if your input sequence file(s) contain assembled data and include the coverage information within each sequence header as described [here](#).

**dereplication** ☐ Remove artificial replicate sequences produced by sequencing artifacts [Gomez-Alvarez, et al, The ISME Journal \(2009\)](#)

**screening**  Remove any host specific species sequences (e.g. plant, human or mouse) using DNA level matching with bowtie [Langmead et al., Genome Biol. 2009, Vol 10, issue 3](#)

**dynamic trimming** ☒ Remove low quality sequences using a modified DynamicTrim [Cox et al., \(BMC Bioinformatics, 2011, Vol. 11, 485\)](#).


Specify the lowest phred score that will be counted as a high-quality base.

Sequences will be trimmed to contain at most this many bases below the above-specified quality.

select

After submission, MG-RAST will process the files. This may take several days, depending upon the number of files submitted, the size of the files, and the number of other files currently waiting in the MG-RAST processing queue.

## SAMSA Pipeline Use: Detailed post-annotation steps



### Step 4: Downloading annotations from MG-FAST

Once MG-FAST has finished processing submitted files, the annotated results can be downloaded using MG-FAST\_API\_downloader.py (part of the SAMSA tools). This tool has the following options:

- Q Quiet mode (no printing to STDOUT); optional
- S Source (RefSeq, UniProt, KEGG, COG, Subsystems, GenBank, etc.); required
- D Data type (Organism, Function, Ontology); required
- A Authorization key, found under "Preferences" at [metagenomics.anl.gov](http://metagenomics.anl.gov); required
- I Annotation ID, found on metagenome page; required

- O Output save file name - should end with .tab; required
- usage Prints usage documentation and exits.

If required flags are not provided when running this tool, it will prompt for the information to be provided.

NOTE REGARDING NEW VERSION OF MG-RAST: To find the Annotation ID, scroll to the bottom of the page; you can click “Sample Data” in the right-hand column menu to jump directly to this portion. The Annotation ID is listed in the table under “Library - metagenome ID”, as shown in the screenshot below:

**MG-RAST**  
metagenomics analysis server

shotgun metagenome

Sam Westrell

Sample Data

category	field	value
env_package	env_package	human-associated
env_package	sample_name	4006_T2_D11_fullFile.merged.cleaned.fastq.gz
library	file_checksum	d6a3b786ed177136e0f9e58f671ee283
library	file_name	4006_T2_D11_fullFile.merged.cleaned.fastq
library	investigation_type	metatranscriptome
library	metagenome_id	4723359.3
library	metagenome_name	4006_T2_D11_fullFile.merged.cleaned.fastq
library	misc_param	Sequences were filtered to remove Illumina PE adaptor contamination.
library	mrna_percent	100
library	rrna_removal_meth	http://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/ribosomal-depletion/ribo-zero/ribo-zero-reference-guide-15066012-01.pdf
library	sample_name	4006_T2_D11_fullFile.merged.cleaned.fastq.gz
library	seq_center	UC Davis DNA Technologies Core, Davis, CA
library	seq_make	HiSeq
library	seq_meth	illumina
library	seq_model	3000

Table of Contents

- Home
- Analysis Statistics
- GSC MixS Info
- DRISSE
- K-mer Profile
- Nucleotide Histogram
- Source Hits Distribution
- Functional Hits
- Taxonomic Hits
- Rank Abundance
- Rarefaction Curve
- Alpha Diversity
- Sequence Length Histogram
- Sequence GC Distribution
- Sample Data**
- Download

Sample MG-RAST\_API\_downloader.py command:

```
$ python MG-RAST_API_downloader.py -S RefSeq -D Organism -I 4628940.3 -A eZYMtTzx5yrYn8pDQx4kUFxC4 -O annotation_results.tab
```

*Result: A downloaded annotation file, containing one line for each submitted, annotated sequence and its best match from the selected database.*

### Step 5: Analyzing annotations

The downloaded annotation file from MG-RAST contains a read-by-read annotation of the submitted metatranscriptome. The next step in the SAMSA pipeline is creating a summary file of this annotation, providing a sorted count of all organisms or transcripts by abundance. This summary file will be imported into R and used for analysis.

To create this summary file for ORGANISM results, use `analysis_counter.py` (part of the SAMSA tools). A sample usage of this script:

```
$ python analysis_counter.py annotation_file -o output_file_name
```

For pipeline use, `-q` will enable quiet mode.

For functional annotations, it is recommended that the MG-RAST internal identifier (M5nr) be preserved. To create this summary file for FUNCTION results, use `analysis_counter.py` with the `"-m"` flag, to preserve M5nr IDs. A sample usage of this script:

```
$ python analysis_counter.py annotation_file -o output_file_name -m
```

*Result: An output file containing a sorted abundance list of all organisms (if using `analysis_counter.py`) or of all functions by M5nr ID (if using the `"-m"` flag). Summary data will be included in the header of the file.*

### **Step 6: Preparing summary files for import into R**

In each output file, the first few lines are summary statistics for the entire file. These lines need to be removed before the file can be imported into R.

To remove these lines, use either:

```
RefSeq_output_reducer.py for an organism file  
func_data_trimmer.py for a functional transcript file
```

For these programs, the input file must be specified with the `"-l"` flag. An output file can be specified using the `"-O"` flag; if this is not included, all files will receive the `"_simplified"` appendix. For more information, run the script with `"-usage"` flag.

*Result: For organism files, `RefSeq_output_reducer.py` will generate a summary file with overall data stripped from the first 6 lines, and with the species names reduced and consolidated down to Genus level. For functional transcript files, `func_data_trimmer.py` will generate a summary file with overall data stripped from the first 6 lines, and with excess unnecessary functional repetition removed.*

### **Optional analysis: Examining different taxonomy levels**

In normal use of this pipeline, organism files are annotated at the Genus level. However, it may be necessary to convert organism results to a higher taxonomic levels for analysis at the Family or Order level. The program `"taxonomy_shifter.py"` is capable of raising the taxonomic level of all entries in a results file (generated in Step 6, above) to a higher taxonomic level.

Also included is a reference database, `Bacteria_Genus_flattened.tsv`, which includes all taxonomic levels for all current bacterial genera. This reference must be specified using the “-R” flag when running the program.

This tool has the following flags needed for usage:

- Q Enables quiet mode
- F Input file, necessary
- R Reference index file, necessary
- T Final taxonomy level desired: Kingdom, Phylum, Class, Order, Family, Genus, necessary
- O Custom output file name (default is `input_file.shifted`)
- V Verbose mode, shows exceptions
- E Exclusion, will exclude all exceptions if present

A sample usage of `taxonomy_shifter.py`:

```
$ python taxonomy_shifter.py -F input_file -R location_of_reference_file -T Phylum -O  
phylum_converted_input_file
```

*Results: A summary file in which all lower taxonomic levels have been raised to the specified taxonomic level and combined into a single result. For example, running the above sample command on an output file will result in a summary output file, as generated in Step 6, but with all entries at the Phylum level.*

## Step 7: Evaluating in R

Once the `analysis_counter` program has run, it will generate output files for each metatranscriptome. These files must be trimmed, using either `RefSeq_output_reducer.py` for organism results, or using `func_data_trimmer.py` for functional results. These files may be converted to a different taxonomic level using `taxonomy_shifter.py`, but are ready for input into R for analysis.

Depending on whether these output files were for organism or transcript reads, they can be analyzed using the R scripts, either `R_organism_script.Rmd` or `R_functional_script.Rmd` (both part of the SAMSA tools).

Both of these R scripts are in Markdown format, with additional comments to explain the analysis being performed by each code block.

Before loading the R analysis programs, ensure that downloaded samples have the following name, either:

- `experimental_filename.output`
- `control_filename.output`

R sorts the imported files into two groups based upon this prefix; if it isn't present, R will not be able to perform pairwise comparison analysis between the two conditional groups.

After files have been imported into R and cleanup steps have been performed (the R program will simplify header names based upon filenames, merge the files into a single table, sort, and simplify naming conventions), further code blocks will allow for the generation either of summary stacked bar graphs, or for the performing of differential analysis between experimental or control conditions.

Further instructions and detailed breakdown of each step in the R analysis script is included in the R markdown scripts.

*Results: The R scripts will generate two stacked bar graphs, showing either the relative or absolute abundance of each different entry across all imported metatranscriptome files. In addition, the DESeq2 section of the R scripts (last code block in the files) will generate a list of all results sorted by adjusted p-value, in tab-delimited form.*

### Further analysis: Isolating functional annotations for a specific organism

When examining a metatranscriptome, researchers may be interested in isolating all functional annotations linked to a specific organism. This step can be performed in the SAMSA pipeline by using “functional\_annotations\_by\_organism.py”, a program that will use the MG-RAST annotation results to screen for all IDs in the functional file that match a specific organism in the organism file.

This tool has the following options needed for usage, specified on the command line:

- N <organism\_name>      The name of the organism you want to search for.
- O <organism\_file>      The RefSeq organism output file.
- F <function\_file>      The RefSeq function output file.
- R <results\_file>      The name given to the results file (outfile).
- Q      Enables quiet mode.
- I <removal\_targets>      A list of targets, such as those below a threshold, to be removed.  
Generated by long\_tail\_threshold.py .
- usage      Prints usage options and exits.

A sample usage of functional\_annotations\_by\_organism.py:

```
$ python functional_annotations_by_organism.py -N Bifidobacterium -O  
control_organism_annotations.tab -F control_function_annotations -R  
Bifidobacterium_functional_summary.tab
```

*Results: functional\_annotations\_by\_organism.py generates a results file with a name based on the “-R” flag. This results file contains a summary of all functional annotations for the organism specified using the “-N” flag. This results file is identical to the file generated after Step 5, although there is no statistics header, so func\_trimmer.py does not need to be run on this file.*

*The results file is in tab-delimited format, and is ready to be imported into R for further analysis (which can be performed using the “R\_functional\_script\_w-md5s.Rmd”).*



### **Further analysis: Converting M5nr internal MG-RAST IDs to RefSeq IDs**

To merge duplicate annotations across multiple databases, MG-RAST uses its own internal identifier system, referred to as M5nr (the “nr” standing for “non-redundant”). Each M5nr internal ID is linked to an external annotation in one or more reference databases used by MG-RAST.

In order to determine the external IDs from other databases that are linked to a specific MG-RAST M5nr ID, an API call can be used to request the relevant information from MG-RAST. For metatranscriptome results, however, it may be necessary to retrieve external IDs for a large number of M5nr IDs at a time. To speed up and automate this process, the program “md5-R-results-query.py” can add external IDs to the results generated by the “R\_functional\_script\_w-md5s.Rmd”.

NOTE: Currently, the M5nr IDs for organisms are not preserved. In order to use this script, the input tab-delimited file must contain the M5nr ID in the right-most column.

Sample usage:

```
$ python md5-R-results-query.py -I input_file_name -O output_file_name
```

Due to the relatively slow speed of retrieving multiple API requests, this program may take several hours to run completely; it can be terminated partly through its run, however, if only the IDs at the beginning of the file are required.

*Results: The output file will appear identical to the input file, saved in tab-delimited form, but will include one extra column on the right, containing the external RefSeq ID corresponding with the MG-RAST M5nr internal ID.*

## Example files and results

SAMSA comes with several example files, designed to test the pipeline and to provide sample results to show how the outputs should appear.

### Sample code-testing data

To see all steps that were performed using these files, please consult the included example bash script, `example_bash_script.bash`. This script breaks down each step, documenting how the commands for each step in the SAMSA pipeline were assembled.

For testing the pre-annotation pipeline, the following files are included in the `sample_files/` folder:

```
control_reads_R1.fastq
control_reads_R2.fastq
experiment_reads_R1.fastq
experiment_reads_R2.fastq
```

(Note that R1 and R2 designate forward and reverse reads, respectively.)

Similarly, for testing the post-annotation pipeline, various files are provided:

After Step 4, annotations downloaded from MG-RAST:

```
control_organism_annotations.tab
control_function_annotations.tab
experiment_organism_annotations.tab
experiment_function_annotations.tab
```

After Step 5, analyzed annotations:

```
control_organism_summary.tab
control_function_summary.tab
experiment_organism_summary.tab
experiment_function_summary.tab
```

After Step 6, summary files prepared for import into R:

```
control_organism_summary_simplified.tab
control_function_summary_simplified.tab
experiment_organism_summary_simplified.tab
experiment_function_summary_simplified.tab
```

Note: because these files are truncated, they aren't suitable for analysis with R. However, the public files can be imported into R and analyzed - see next page.

## **Playing with real-life public metatranscriptome data**

Ah, but what about some real life data? How does SAMSA handle full metatranscriptomes?

We provide two public metatranscriptomes, designated as 4012 and 4015. The raw data from these metatranscriptomes is available on NCBI's Short Read Archive (SRA), here:

<http://trace.ncbi.nlm.nih.gov/Traces/sra/sra.cgi?study=SRP071017>

The MG-RAST RefSeq annotations of these files were downloaded, both for organism and for functional annotations, and analyzed using the `analysis_counter.py` program. Their outputs are provided in the `public_output_files/` folder:

`4012_P1_RefSeq_organism.tab.output`

`4012_P1_RefSeq_function.tab.output`

`4015_P1_RefSeq_organism.tab.output`

`4015_P1_RefSeq_function.tab.output`

These output files were then prepared for input into R, using the `RefSeq_output_reducer.py` program on the organism files, and the `func_data_trimmer.py` program on the functional annotation files, to create the following files:

`4012_P1_RefSeq_organism.tab.output_simplified`

`4012_P1_RefSeq_function.tab.trimmed.output`


`4015_P1_RefSeq_organism.tab.output_simplified`

`4015_P1_RefSeq_function.tab.trimmed.output`


These files are split into two output folders, labeled “organism” and “function”, in `/public_output_files/public_files_ready_for_R/`, and are ready to be imported into R for analysis.

Results of data file analysis, from public data, as exported by R


Organism results, relative abundance graph (as a percentage of total sample):




Organism results, total abundance graph (by raw annotation counts):



Functional results, relative graph (as a percentage of total sample):



Functional results, total abundance graph (by raw annotation counts):



So, looking at these graphs, we can see that both individuals have relatively similar metatranscriptomes in terms of organism abundance. However, these organisms are expressing different pathways in the different individuals, as we see in the functional annotations.

Note that R can also run DESeq analysis on these datasets, although that is not included in the sample data.

## Frequently asked questions (FAQs)

Questions will be added here as needed.

### 1. How do I download all of the pipeline scripts and files from Github?

There are two different methods for downloading the SAMSA pipeline and associated files. From the Github page (<https://github.com/transcript/SAMSA>), click the green button in the upper right, labeled “Clone or download”. From here, you may either download the entire pipeline in a zipped folder, or you can use Git’s clone function to copy the pipeline to your computer.

If you choose to use the clone function, copy the clone URL, and then open your terminal. Use the following command:

```
# git clone https://github.com/transcript/SAMSA.git
```

This will create a folder containing all SAMSA programs and files at the home location of your computer (~/ if using a Mac, within the Cygwin environment if using Windows).

### 2. Where do I generate my MG-RAST authorization key?

After creating an account on MG-RAST, the authorization key can be generated under Account Management -> Preferences. When logged in, select the small green person icon at the top right of the page, and then click to manage your personal preferences. Under “Web Services”, click the button labeled “generate new key”. Note that you must also click “set preferences” at the bottom of the screen to save this newly generated key.

Note that a generated webkey is only valid for a limited time, usually 7 days. After this time, you must generate a new webkey to use in your commands.

### 3. How do I generate log files to preserve my steps and commands?

Generating log files with each step in the SAMSA pipeline is straightforward. Normally, SAMSA prints the commands used, as well as details about the processing steps, to STDOUT (unless quiet mode is enabled). These details and information can be redirected to a log file by amending the command used. For example, for downloading MG-RAST annotations:

```
# python SAMSA/MG-RAST_API_downloader.py -S RefSeq -A @@AuthorizationKey@@ -D Organism -I @@AnnotationID@@ -O control_organism_annotations.tab > api_download.log
```

In this example, the messages and announcements will be saved to a log file named “api\_download.log”.

#### **4. What version of Phred scores should I use?**

Phred quality scores are measures of the potential error rate for calling each nucleotide base when sequencing a read. Sanger sequencing and Illumina post-1.8 format sequencing produces Phred scores in +33 notation, while pre-1.8 format Illumina sequencing produces Phred scores in +64 notation.

Both of these cases are acceptable for the SAMSA pipeline; the Trimmomatic read cleaning step automatically detects whether the FASTQ input file is in Phred+33 or Phred+64 notation, and adjusts accordingly. Most likely, the Phred scoring will be in +33 notation, as this is the current scoring notation used by Sanger, Illumina, and for storage in the NCBI archive.



## Contact information

Every piece of software is a constant work in progress; SAMSA is no different. There are likely numerous bugs hidden in the code. If you encounter an issue and require assistance, you can contact me at:

Sam Westreich

[stwestreich@ucdavis.edu](mailto:stwestreich@ucdavis.edu)

In addition, as I encounter frequently asked questions, I will add the answers into the documentation.